


Bellacide® 355 not just a biocide

BETTER MICROBIOLOGICAL PERFORMANCE

Bellacide 355 as a Azole substitute

Reduce the azoles in your formulations by 50%, while improving your copper corrosion rates




Equipment: Rotating coupon test baths
 Coupons: Copper, grade Cu 106 (99.9% min Cu)
 Duration: Rotating at 150 rpm for 90 days
 Bath temperature: 30°C (86°F)

Water make up
 60ppm Ca²⁺
 18ppm Mg²⁺
 41.4ppm Cl⁻
 427ppm HCO₃⁻
 200ppm SO₄²⁻

Bellacide 355 as a Corrosion inhibitor

Corrosion rate of N80 steel in mpy after 72 hours



Our wide selection of chemistry can be tailored for any treatment scenario.

Learn more at www.wateradditives.com